mindray

M09G

Tympanic Thermometer

Instructions for Use

English 1
Deutsch 17
Español 35
Francais 53
Italiano 71
Nederlands 89
Portugês 107
Svenska 125
Türkçe 141

046-023413-00(5.0)

Statement

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. (hereinafter called Mindray) owns the intellectual property rights to this product and this manual. Disclosure of the information in this manual in any manner whatsoever without the written permission of Mindray is strictly forbidden.

This manual provides the instructions necessary to operate the product in accordance with its function and intended use. Observance of this manual is a prerequisite for proper performance and correct operation, and ensures patient and operator safety.

Mindray is responsible for safety, reliability and performance of this product only in the condition that:

- This product is operated under strict observance of this manual.
- This product is not damaged by human factors. Human factors refer to unintentional falling, intentional damaging, etc.

In the event that it becomes necessary to return a unit to Mindray, please contact the Mindray Service Department and obtain a Mindray Customer Service Authorization Number. The Mindray Customer Service Authorization Number must appear on the outside of the shipping container. Return shipments will not be accepted if the Mindray Customer Service Authorization Number is not clearly visible. Please provide the model number, serial number, and a brief description of the reason for return. The customer is responsible for freight charges when this product is shipped to Mindray for service (including any relevant customs fees or other freight related charges).

mindray, MINDRAY, MINDRAY and are the registered trademarks or trademarks owned by Mindray in China.
© 2022-2026 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved.

- Release time: 2026-06
- Revision: 5.0

WARNING

- If there is any inconsistency or ambiguity between the English version and this version, the English version shall prevail.

Company Contact

Manufacturer:	Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Address:	Mindray Building, Keji 12th Road South, High-tech Industrial Park, Nanshan, Shenzhen 518057, P.R.China
Website:	www.mindray.com
E-mail:	service@mindray.com
Tel:	+86 755 81888998
Fax:	+86 755 26582680
EU-Representative:	Shanghai International Holding Corp. GmbH (Europe)
Address:	Eiffestraβe 80, 20537 Hamburg, Germany
Tel:	0049-40-2513175
Fax:	0049-40-255726

Notification of Adverse Events

As a health care provider, you may report the occurrence of certain events to SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD., and possibly to the competent authority of the Member state in which the user and / or patient is established.

These events, include device-related death and serious injury or illness. In addition, as part of our Quality Assurance Program, SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. requests to be notified of device failures or malfunctions. This information is required to ensure that SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. provides only the highest quality products.

Safety

Safety Information

WARNING

- Check the thermometer, cradle and cable before using. Do not use them for temperature measurement in case of damage, block, deterioration or visible contamination.
- To avoid explosion hazard, do not use the equipment in the presence of oxygen-rich atmospheres, flammable anesthetics, or other flammable agents.
- The equipment is not intended to be used within the Magnetic Resonance (MR) environment.
- Magnetic and electrical fields are capable of interfering with the proper performance of the equipment. For this reason make sure that all external devices operated in the vicinity of the equipment comply with the relevant EMC requirements. Mobile phones, X-ray equipment or MRI devices are possible sources of interference as they may emit higher levels of electromagnetic radiation.
- To reduce the influence of high-frequency surgical procedure on the measurement, keep the thermometer away from the cables of the electrosurgery unit (ESU).
- The thermometer may not meet the performance specifications if stored or used outside the specified temperature and humidity ranges. If the performance of the device is degraded due to aging or environmental conditions, contact your service personnel.
- Always install or carry the equipment properly to avoid damage caused by drop, impact, strong vibration or other mechanical force.
- The equipment should be installed by authorized Mindray personnel.
- Connect only approved devices to this equipment. Devices connected to the equipment must meet the requirements of the applicable IEC standards (e.g. IEC 60950 safety standards for information technology equipment and IEC 60601-1 safety standards for medical electrical equipment). The system configuration must meet the requirements of the IEC 60601-1 medical electrical systems standard. Any personnel who connect devices to the equipment's signal input/output port are responsible for providing evidence that the safety certification of the devices has been performed in accordance to the IEC 60601-1. If you have any questions, contact Mindray.
- Use only accessories specified in this manual. Using any accessory other than the stated in the accessory section would lead to failed or inaccurate measurements.
- Use disposable probe covers for temperature measurement. Failure to use a probe cover can cause inaccurate temperature readings, and patient cross-contamination.
- The probe cover is intended for use with only Mindray specified equipment.
- The probe cover is disposable. Re-use of probe cover may result in patient cross-contamination.
- Check the disposable probe cover for damage before using. Never use any probe cover for temperature measurement in case of damage or contamination.
- Check the compatibility between the equipment and accessory before use. Otherwise, product performance may be compromised or personal injury may occur.
- This equipment is used for single patient at a time.
- The equipment is not water-proof. Keep it away from water.
- Do not apply extra force on the patient's ear when performing measurements. Damage to the ear canal could result.

- Do not take measurements on bleeding or damaged skin.
- Do not take measurements on a skin where inflammation occurs.
- Do not take measurements if there is bleeding or pus in the ear canal.
- Do not take measurements on a patient who is under in-ear treatment.
- If a patient has a deformation of the ear canal or face, measurements may fail.
- Wash hands after temperature is taken. This will significantly reduce the risk of cross contamination and nosocomial contamination.
- The measurements provided by the thermometer should not be used as the only basis for diagnosis or therapy decisions. They must be used in conjunction with clinical signs and symptoms. Misinterpreting measured values or other parameters may result in patient hazards.
- At the end of its service life, the equipment, as well as its accessories, must be disposed of in compliance with the guidelines regulating the disposal of such products.
- For disposal of accessories, follow local regulations or your hospital's policy regarding disposal of such accessories. Do not dispose randomly.

NOTE

- For TrueTymp cradle independent, the equipment uses the mains plug of the adapter as isolation means to the mains power. Do not locate the equipment or adapter in a place difficult to operate the mains plug.
- The operator should stand in front of the device during normal use.
- The equipment software copyright is solely owned by Mindray. No organization or individual shall resort to modifying, copying, or exchanging it or to any other infringement on it in any form or by any means without due permission.

Equipment Symbols

Symbol	Description	Symbol	Description
	Do not reuse		Caution
	Manufacturer		Date of manufacture
LOT	Batch code		Use by date
	Type BF applied part	SN	Serial number
	Temperature limit		Refer to instruction manual/booklet (Blue)
	Atmospheric pressure limita tions		Humidity limitations
	Take measurement		Timer
	Direct current		Alternating current
	Non-ionizing electromagnetic radiation		Recycle
	Class II Equipment	UDI	Unique device identifier
MD	Medical Device	PA	Pattern Approval Certification
EU REP	Authorized representative in the European Union		Dispose of in accordance to your country's requirements

Symbol	Description	Symbol	Description
Rx ONLY	Federal law (USA) restricts this device to sale by or on the order of a physician.		
	Complies to Administrative Measure on the Control of Pollution Caused by Electronic Information Products in China. Contains substances with 10 years environmentally friendly use period (EFUP).		
TrueTymp™	TrueTymp™ is Mindray's new ear temperature measurement algorithm using infrared energy from the tympanic membrane, which can measure body temperature accurately, quickly and conveniently.		
	This product is provided with a CE marking in accordance with the regulations stated in Regulation (EU) 2017/745 of the European parliament concerning Medical Devices. The number adjacent to the CE marking (0123) is the number of the EU-notified body certified for meeting the requirements of the Regulation. Note: The product complies with the Council Directive 2011/65/EU.		

Intended Use

Intended Purpose

The Tympanic Thermometer is intended for spot checking, displaying, reviewing, storing and transferring temperature of single patient.

NOTE

- According to the conclusion of clinical evaluation and residual risk evaluation, for the intended patients, there is no known side effects that can occur during or after the use of the medical device. And there is no need for the operator to make extra preparations. Thus, no residual risk associated with using the medical device should be disclosed due to the risk management report.

Indication for Use

The Tympanic Thermometer is intended to be used for spot checking the temperature on adult, pediatric, and neonatal patients in healthcare facilities by clinical physicians or appropriate medical staff under the direction of physicians.

Intended Users

The Tympanic Thermometer is intended for use only by clinical professionals or under their guidance. It must only be used by persons who have received adequate training in its use. Anyone unauthorized or untrained must not perform any operation on it.

Intended Patient Population

The Tympanic Thermometer is applicable for adult, pediatric, and neonatal patients.

Intended Medical Conditions

The Tympanic Thermometer is expected to be used in medical institutions, and its application fields include: emergency care, internal medicine and surgical care, neonatal care, elderly care, Ambulance Surgical Center, postoperative recovery, Physician Office.

Limitations

N/A.

Contra-indications

Inflammation, trauma, postoperative and other local lesions at the measurement site.

Side-effects

None.

Through clinical data from literature and clinical data from post-market surveillance activity of subject device, there is no side-effects identified. After search the literature of similar device, the results of SOTA evaluation shown that undesirable effects may include myocardial damage.

Unpacking and Checking

Before unpacking, examine the packing case carefully for signs of damage. If any damage is detected, contact the carrier or us.
If the packing case is intact, open the package and remove the equipment and accessories carefully. Check all materials against the packing list and check for any mechanical damage. Contact us in case of any problem.

Power on the Thermometer

Before first-time use, you may need to power on the thermometer as instructed below:

1. Connect the cradle with correct power source.
 - For TrueTymp cradle independent, connect the cradle with an AC adapter via the USB connector. Then connect the AC adapter to AC power.
 - For TrueTymp cradle with link, turn on the patient monitor.
2. Place the thermometer to the cradle.

If the screen on the thermometer lights up, the thermometer is ready for use.

Product Overview

The Tympanic Thermometer consists of the ear thermometer, cradle (including cradle independent and cradle with link), probe cover, power adapter and battery.

Tympanic Thermometer

1. Eject button: press the button to eject the disposable probe cover from the probe tip
2. NFC (Near Field Communication) tag: used to transmit the temperature measurements
3. Communication and power input contacts
4. Security tether
5. Measurement key: press the key to start a temperature measurement.
6. Screen: displays the measurements and system information, including battery status, measurement unit etc.
7. Unit key: press and hold the key to convert the measurement unit to °C (Celsius) or °F (Fahrenheit).
8. Timer key:
 - Press and hold the key to start a 60s timer
 - In the process of timing, press the key to quit
9. Probe: Applied part
10. Review key:
 - Press the key to enter Review screen and view the latest measurement
 - Under Review screen, press the key to show earlier measurements
 - Under Review screen, press and hold the key to exit Review mode

TrueTymp Cradle

TrueTymp cradle independent TrueTymp cradle with link

1. Probe cover box
2. Probe covers
3. USB connector: used to connect an AC adapter
4. Probe cover box support
5. Communication and power connector: used to connect a Mindray monitor

On-screen Symbols

	Battery icon: displays in the top left corner to indicate the current charge of the battery. When the solid parts fills the symbol dynamically, the battery is charging.
	Battery icon: displays in the top left corner. When it flashes, the battery is low or nearly depleted. Charge the battery immediately.
	Install probe cover icon: indicates that a probe cover needs to be installed before taking measurements.
	Eject probe cover icon: indicates that the measurement is finished, and the probe cover can be ejected.

Working Modes


The thermometer supports the following two working modes:

- Adjusted mode: in this mode, the thermometer measures the tympanic temperature, and, after compensation, displays an equivalent oral temperature. Adjusted mode is the intended working mode for all users.
- Direct mode: the thermometer measures and displays the measured temperature. Direct mode is available only for authorized service personnel.

Preparing the Probe Covers

For easy access to probe covers when taking measurements, you need to install a box of probe covers onto the cradle. Follow this procedure:

1. Tear the probe cover box open as indicated.
2. With the opening facing up, position the box onto the probe cover box support in the cradle.

3. Place the thermometer back to the cradle.

Factors Affecting Measurement Accuracy

The measurements may be inaccurate or impossible in the following situations:

- Disposable covers other than those specified in this manual are used.
- The disposable cover used is damaged, cracked or dirty.
- The probe was not pointed right at the ear canal or shaken during measurement.

- The probe lens window is dirty or the probe is damaged.
- The thermometer is used or stored outside the specified temperature and humidity ranges.
- The patient is wearing a hearing aid or earbud.
- The patient's ear has been squeezed because of side-lying.
- The patient is in an agitated emotional state.
- There is excessive patient movement.
- There are excessive earwax or oil in the ear canal.

Performing Measurement

Before measurement, make sure all warnings listed in the previous section have been understood.

WARNING

- Follow the steps below strictly when taking measurements. Otherwise, the measurement may fail or be inaccurate.


NOTE

- If the thermometer is moved from one place to another of a great temperature difference, wait for at least 30 minutes before performing measurement on a patient.
- Allow a patient from a hot or cold environment to adjust to the current room temperature before taking measurements.
- Wait for at least an one minute interval between performing measurements on the same site of a patient.


1. Take the thermometer out from the cradle. Then the thermometer starts a selftest. During this process, observe the screen, and make sure the numerics are displayed as follows. If not, remove the thermometer from use.


 If the test is passed, the screen displays as follows, and the thermometer is ready for measurement.


2. Press the tip of the probe into a probe cover from the container until you hear a click.

3. Lift the auricle of the patient up to the back with one hand, and place the tip of the probe into the ear canal of the patient with the other hand. Make the sensor part stick closely to the ear canal and the probe window point directly at the ear drum. Then press the measurement key.


4. Wait till you hear a beep, which indicates the completion of the measurement. Check the screen for the temperature reading.
5. When the measurement is finished, point the tip of the probe to an specialized waste receptacle, and press the eject button to eject the probe cover.
6. Upload measurements and place the thermometer back to the cradle.
 - For thermometers used with a TrueTymp cradle Independent: put the NFC tag close to the NFC device to upload the measurement until you hear a beep. Then place the thermometer back to the cradle.
 - For thermometers used with a TrueTymp cradle with link: place the thermometer back to the cradle and the measurement is automatically updated on the monitor where the cradle is installed.

7. Check the monitor screen and make sure the measurement is updated.

Reviewing History Measurements

You can press the review key to check the history measurements. Up to 20 measurements can be viewed.

1. Temperature reading
2. Temperature unit
3. Number of the current measurement

The measurement displays for 10 seconds if no operation is detected. To view the measurements, press the review key.


Exiting the Standby Mode

After a measurement completes, if no operation is detected in 10 seconds, the thermometer enters standby mode automatically.

Press any key to exit the standby mode.

Using the Timer

The thermometer provides a 60 second counting-up timer. Press and hold the timer key to start timing. A beep is given every 15 seconds and the timer stops once the 60th second is reached. During the timing, you can also press the timer key to stop the timer.

1. Counting seconds
2. Progress bar

Changing Temperature Unit

On the thermometer, temperature can be displayed in either °C or °F. Press and hold the unit key to switch the temperature unit from the current one to the other.

Cleaning and Disinfection

WARNING

- Never immerse any part of the thermometer or cradle in liquids or allow liquid to enter the interior.
- Do not clean or disinfect the metal parts of this product. Otherwise, metal part corrosion or product reliability deterioration could result.
- Dirt on the probe window may affect the measurement accuracy. If needed, clean the window with a medical cotton bud moistened with 70% ethanol. Remove the liquid residue on the window with a clean cotton bud if needed.
- We recommend that the product be disinfected only when necessary as determined by your hospital's policy, to avoid long term damage to the product.
- Use only cleaners, disinfectants and methods specified in this manual. Using unapproved substances or methods may damage the equipment and void the warranty.

Recommended Cleaning and Disinfecting Agents

For the recommended disinfectants, refer to the Appendix.

Cleaning

1. Clean the thermometer and cradle with cotton or soft cloth moistened with water.
2. After cleaning, wipe off the water with a soft cloth.
3. Allow the thermometer and cradle to air dry.

Disinfection

1. Clean the thermometer and cradle as instructed above.
2. Disinfect the thermometer and cradle with cotton or soft cloth moistened with one of the recommended disinfectants.
3. After disinfection, be sure to wipe off the disinfectant left on the thermometer and cradle with a soft cloth moistened with water.
4. Allow the thermometer and cradle to air dry.

Maintenance

WARNING

- This product contains no user serviceable parts. No modification is allowed. Contact your service personnel in case of any problem.
- Do not open the equipment housings. All servicing and future upgrades must be carried out by the personnel trained and authorized by our company only.
- The service personnel must be properly qualified and thoroughly familiar with the operation of the equipment.

NOTE

- If needed, contact the manufacturer for circuit diagrams, component part lists, descriptions, calibration instructions, or

other information concerning the repair of the product.

Battery

The thermometer is powered by an embedded rechargeable lithium battery. The battery is automatically charged if placed onto a cradle that has been connected to power supply.

WARNING

- Keep the battery charge observed and charge the battery timely.
- Use only the specified adapter for power supply and charging. Otherwise, equipment damage or patient injury may result.
- If the thermometer is to be stored for an extended long time, before storing, charge the battery to at least two bars. Besides, charge the battery every three months during the storage. Failing to do so could damage the thermometer.
- The lithium-ion battery has a service life of five years. Contact your service personnel if the battery reaches the end of its service life or the performance deteriorates greatly.

Verification

Verification should be performed if you suspect that the measurement values are incorrect. Besides, the thermometer should be verified for at least once a year. Contact the service personnel for verification.

Troubleshooting

This section lists the problems that are likely to occur and possible solutions. If the problem persists after corrective actions have been taken, contact your service personnel.

Symptom	Possible cause	Attemptable solution
The screen cannot be woken up	The battery charge is depleted. The thermometer failed.	Charge the battery immediately. Contact the service personnel.
	The ambient temperature is too low, exceeding applicable operating range.	Move the thermometer and patient to an environment that meets the specified temperature and humidity ranges. Wait for 30 minutes and take measurement again.
	The ambient temperature is too high, exceeding applicable operating range.	
	The temperature exceeds the lower limit of the measurement range.	Replace the probe cover and take measurement again. The patient's temperature may be too high or too low. Check the patient condition.
	The temperature exceeds the upper limit of the measurement range.	
	Patient motion. The thermometer failed.	Ask the patient to relax and stay still. Take measurement again. The selftest failed. Contact the service personnel.

Accessories

The accessory material that contacts the patients or other staff has undertaken the biocompatibility test and is verified to be in compliance with ISO 10993-1.

WARNING

- Use accessories specified in this manual. Using other accessories may cause damage to the equipment or not meet the claimed specifications.
- Single-use accessories are not designed to be reused. Reuse may cause a risk of contamination and affect the measurement accuracy.

CAUTION

- Check the accessories and their packages for any sign of

damage. Do not use them if any damage is detected.

- Use the accessories before the expiry date if their expiry date is indicated.
- The disposable accessories shall be disposed of according to hospital's regulations.

Probe cover

Model	Usage	Patient category	Description	Part No.
ET-001	Disposable	Adult, pediatric, neonate	Cover, 200 pcs/pack	040-006890-00

Others

Part No.	Description
115-085932-00	TrueTymp cradle independent, MR495
115-085931-00	TrueTymp cradle with link, MR496
115-089747-00	AC adapter, power cord (Europe), MANGO18SC-05
115-089746-00	AC adapter, power cord (UK)
115-089748-00	AC adapter, power cord (U.S.A.)
115-089749-00	AC adapter, power cord (Korea)
115-089750-00	AC adapter, power cord (Australia)
115-089751-00	AC adapter, power cord (Brazil)
115-089752-00	AC adapter, power cord (India)
022-000644-00	AC adapter, power cord (Switzerland)
040-007149-00	Security tether

Specifications

Classifications

The equipment is classified, according to IEC60601-1:

Type of protection against electrical shock	Class II
Degree of protection against electrical shock	Type BF applied part
Degree of safety of application in the presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide	The equipment is not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide
Defibrillation-proof applied part	None
Mode of operation	Continuous

Environmental Specifications

Item	Temperature	Relative humidity	Barometric
Operating conditions	10°C to 40°C	10% to 95%, noncondensing	427.5 mmHg to 805.5 mmHg (57.0 kPa to 107.4 kPa)
Storage and transport conditions	-25°C to 60°C	10% to 95%, noncondensing	120 mmHg to 805.5 mmHg (16.0 kPa to 107.4 kPa)

Power Supply Specifications

Battery (embedded in the thermometer)

Battery type	Rechargeable lithium-Ion battery, LI11S002G, 3.6V, 3200mAh
Run time	At least 4500 times of measurement when the thermometer is powered by a new fully-charged battery at 25 °C±5 °C (Process of a typical measurement: take the thermometer from the cradle, take a measurement, and upon completion, put the thermometer back to the cradle)

TrueTymp cradle independent (with AC adapter)

Model	MANGO18SC-05
AC power input	100 to 240 VAC (±10%), 50/60Hz (±3Hz), 0.5A Max

TrueTymp cradle with link

DC power input	5 VDC±5%, 0.5A, powered by the monitor

Physical Specifications

Item	Model	Weight	Depth	Width	Height
Tympanic thermometer	M09G	≤160g	151±2mm	36±2mm	39±2mm
TrueTymp cradle independent	MR495	≤110g	166±2mm	42±2mm	57±2mm
TrueTymp cradle with link	MR496	≤110g	166±2mm	42±2mm	57±2mm

Hardware Specifications

Screen type	Segment LCD
Screen size	1.3 inch

Interface Specifications

Tympanic thermometer	1 power input and communication interface 1 NFC interface
TrueTymp cradle independent	1 power interface 1 Type C USB connector
TrueTymp cradle with link	2 power input and communication interfaces

NFC Specifications

Protocol	ISO/IEC14443A
Operating frequency	13.56 MHz

Measurement Specifications

Standards	Meet the standard of ISO 80601-2-56	
Display range	20.0°C to 42.2°C	
Measurement range	20.0°C to 42.2°C	
Accuracy	With probe cover	Without probe cover
	At a room temperature of 16°C to 40°C: • ±0.2°C in measurement range 35.0°C to 42.2°C; • ±0.3°C in measurement range 20.0°C to 35.0°C (not included). At a room temperature of 10°C to 16°C (not included): ±0.3°C	At a room temperature of 16°C to 40°C: • ±0.1°C in measurement range 35.0°C to 42.2°C; • ±0.2°C in measurement range 20.0°C to 35.0°C (not included). At a room temperature of 10°C to 16°C (not included): ±0.2°C
Resolution	0.1°C	
Measurement time	≤2s	

Statistical Results of Clinical Investigation Data

Measurement site	Reference body site	Clinical BIAS (Δcb)	Limits of Agreement (LA)	Clinical Repeatability (σr)
Ear	Oral	-0.03	0.31	0.10

Operating Environment

Operating system	FreeRTOS 9.0.0
CPU	STM32F411RCTx
RAM	128 KB
Network	NFC, RS232

EMC and Radio Regulatory Compliance

The device meets the requirements of IEC 60601-1-2: 2014.

NOTE

- Use of this device adjacent to or stacked with other device should be avoided because it could result in improper operation. If such use is necessary, this device and the other device should be observed to verify that they are operating normally.
- The non-ME EQUIPMENT (e.g. ITE) that is a part of an ME SYSTEM may be disrupted by the electromagnetic interference of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the non-ME EQUIPMENT or shielding the location.
- Use of this device adjacent to or stacked with other device should be avoided because it could result in improper operation. If such use is necessary, this device and the other device should be observed to verify that they are operating normally.
- This device is intended for use in professional healthcare facility EMC environment. If it is used in special environment, such as magnetic resonance imaging environment, the equipment/ system may be disrupted by the operation of nearby equipment.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of this device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this device could result.

Guidance and Declaration – Electromagnetic Emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
Emission tests	**Compliance**	**Electromagnetic environment – guidance**
Conducted and radiated RF EMISSIONS CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic device.
Conducted and radiated RF EMISSIONS CISPR 11	Class A	The device is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes
Conducted and radiated RF EMISSIONS CISPR 11	Class B	The product is suitable for use in all establishments including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic distortion EMISSIONS IEC61000-3-2	Class A	
Voltage Fluctuations/ Flicker EMISSIONS IEC 61000-3-3		

If the equipment is operated within the electromagnetic environment listed in Table Guidance and Declaration - Electromagnetic Immunity, the system will remain safe and provide the following essential performance:

- Operating mode
- Accuracy
- Function

Guidance and Declaration - Electromagnetic Immunity
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge(ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	±1 kV line(s) to line(s) ±2 kV line(s) to earth	
Voltage dips and Voltage interruptions IEC 61000-4-11	0% U_T for 0.5 cycle 0% U_T for 1 cycle and 70% U_T for 25/ 30 cycles 0% U_T for 250/300 cycles	0% U_T for 0.5 cycle 0% U_T for 1 cycle and 70% U_T for 25/ 30 cycles 0% U_T for 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of our product requires continued operation during power mains interruptions, it is recommended that our product be powered from an uninterruptible power supply or a battery.
RATED power frequency magnetic fields IEC 61000-4-8	30 A/m 50 Hz/60 Hz	30 A/m 50 Hz/60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note: U_T is the AC mains voltage prior to application of the test level.			

Guidance and Declaration – Electromagnetic Immunity			
The device is intended for use in the specified electromagnetic environment. The customer or the user of the device should assure that it is used in such an environment as described below.			
Immunity test	**IEC60601 test level**	**Compliance level**	**Electromagnetic environment – guidance**
Conducted disturbances induced by RF fields IEC61000-4-6	10 Vrms 150 kHz to 80 MHz	10 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2\sqrt{P}$ 80 MHz ~ 800 MHz $d = 1.2\sqrt{P}$ 800MHz ~ 2.5GHz $d = 2.3\sqrt{P}$
	6 Vrms in ISM bands and amateur radio bandsa between 0,15 MHz and 80 MHz	6 Vrms	

Radiated RF EM fields IEC61000-4-3	3V/m 80 MHz to 2.7 GHz	3V/m	where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey[a], should be less than the compliance level in each frequency range[b]. Interference may occur in the vicinity of equipment marked with the following symbol: 
Proximity fields from RF wireless communic ations equipment IEC61000- 4-3	27 V/m 380–390 MHz	27 V/m	
	28 V/m 430–470 MHz, 800–960 MHz, 1700–1990 MHz, 2400–2570 MHz	28 V/m	
	9 V/m 704–787 MHz, 5100–5800 MHz	9 V/m	

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

[a]. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.
[b]. Over the frequency ranges 150 kHz to 80 MHz, field strengths should be less than 3V/m.

Recommended Separation Distances between Portable and Mobile RF, Communications Equipment and This Equipment

The equipment is intended for use in an electromagnetic environment in which radiated RF disturbance are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communication equipment.

Rated Maximum Output power of Transmitter Watts (W)	**Separation Distance According to Frequency of Transmitter (m)**		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters at a maximum output power not listed above, the recommended separation distanced in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.
Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Radio Regulatory Compliance

The radio device used in this product is in compliance with the essential requirements and other relevant provisions of Directive 2014/53/EU.

This page intentionally left blank.

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Appendix: List of Recommended Disinfectants

The recommended disinfectants are as follows:

Product name	Product type	Manufacturer
Water	Liquid	/
Ethanol, 70%	Liquid	/
N-propanol (1-propanol), 50%	Liquid	/
Alpet® D2 Surface Saniyizing Wipes	Wipes	BEST SANITIZERS INC™.
Clorox Dispatch® Hospital Cleaner Disinfectant Towels with Bleach	Wipes	Clorox professional products company
Clorox Healthcare® Bleach Germicidal Wipes	Wipes	
Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipe	Wipes	
Diversey Oxivir® TB Wipes	Wipes	Diversey Inc
Hydrogen peroxide, 3%	Liquid	/
Isopropanol, 70%	Liquid	/
Metrex CaviCide1™	Liquid, spray	METERX® RESEARCH
Metrex CaviWipes™	Wipes	
PDI Sani-Cloth® AF3 Germicidal Disposable Wipe	Wipes	PDI Inc.
PDI Sani-Cloth® Bleach Germicidal Disposable Wipe	Wipes	
PDI Sani-Cloth® Plus Germicidal Disposable Cloth	Wipes	
PDI Super Sani-Cloth® Germicidal Disposable Wipe	Wipes	
Perform® Classic Concentrate OXY, 0.5%	Powder	Schülke & Mayr GmbH
Rely+On™ Virkon® High Level surface Disinfectant, 1%	Powder	Antec International Ltd
Sodium hypochlorite bleach 0.5%	Liquid	/
Virex® II 256 (1:256)	Liquid, spray	Diversey Inc
Virex® TB	Liquid	
Health Essence Disinfecting Effervescent Tablets	Tablets	Beijing ChangJiangMai Medical Science Technology Co. Ltd
Environmental surface of Health Essence brand Disinfectant	Liquid, spray	
Health Essence Double chain quaternary ammonium salt Disinfectant	Liquid	
DIAN'ERKANG® surface wipes	Wipes	Shanghai Likang Disinfectant Hi-Tech Co., Ltd
DIAN'ERKANG® Surface Disinfectant	Liquid	
DIAN'ERKANG® Disinfectant Spray	Liquid, spray	
Clinell Universal Wipes	Wipes	GAMA Healthcare Ltd
Clinell® Sporicidal Wipes	Wipes	
Tristel Duo™	Liquid, spray	Tristel solutions Limited
Tristel Jet	Liquid, spray	
Tristel Fuse For Surfaces, 196ppm	Liquid	
Surfanios Premium, 0.25%	Liquid	ANIOS LABORATORIES
Anios Surfa 'safe	Liquid, spray	
Wip' Anios EXCEL	Wipes	
Aniosurf ND premium, 0.25%	Liquid	
Mikrobac® Tissues	Wipes	BODE Chemie GmbH
Cleanisept® Wipes	Wipes	Dr. Schumacher GmbH

Product name	Product type	Manufacturer
Mikrozid® PAA Wipes	Wipes	Schülke & Mayr GmbH
Mikrozid® Sensitive Wipes	Wipes	
Glutaraldehyde, 2%	Liquid	/
Ecolab Incidin® OxyWipe S.	Wipes	Ecolab Deutschland GmbH
Descosept® forte	Liquid	Dr. Schumacher GmbH
Descosept® AF	Liquid	
Dismozon® plus, 0.4%	Powder	BODE Chemie GmbH
mikrozid® AF Wipes	Wipes	Schülke & Mayr GmbH
Terralin® Liquid	Liquid	

CE 0123

EU REP

Shanghai International Holding Corp. GmbH(Europe)
Eiffestraβe 80, 20537 Hamburg Germany

www.mindray.com
Mindray Building, Keji 12th Road South, High-tech
Industrial Park, Nanshan, Shenzhen 518057, P.R. China
Tel: +86 755 81888998
Fax:+86 755 26582680
E-mail: service@mindray.com